



THE SPECIALIST
MANUFACTURER
OF HIGH QUALITY
ANAESTHETIC
VAPORISERS

MSS
INTERNATIONAL LTD

PRECISION
ENGINEERING
AND SERVICING
WHERE IT MATTERS



Founded in 1996, MSS International Limited has built a successful company, both in the manufacturing, developing and servicing of anaesthetic vaporisers.

Directors John and Moira McInnes have over 60 years combined experience in the industry, making the name of MSS synonymous with quality products and service.

MSS developed their first Tec 3 anaesthetic vaporiser in April 1997 and since then have manufactured nearly 20,000 units, exporting to markets in the US, Europe, China, Africa and the Middle East. We are truly a global business.

Today the Tec 3 has seen many developments and improvements and the latest MSS3 is engineered to the highest standards, with 90% of parts now manufactured and all vaporisers are assembled on site in its own dedicated machine and workshop facility.

**DELIVERING
HIGH STANDARDS
LOCALLY AND
INTERNATIONALLY...**

At MSS attention to detail is important. Working with top grade materials, worked on by specialist engineers and skilled staff.

The resulting product meets and exceeds industry regulations and criteria. Making MSS one of the market leaders in anaesthetic vaporiser manufacture, with MSS vaporisers in operation all over the world.

MSS are BSI certified and CE accredited, supporting the dedication to provide a product which delivers high standards locally and internationally.



THE SPECIALIST MANUFACTURER OF HIGH QUALITY ANAESTHETIC VAPORISERS

MSS
INTERNATIONAL LTD

THE LATEST VAPORISER TECHNOLOGY...

Over the years MSS has developed and improved the MSS3, creating a stable vaporiser with output precision and a longer duration of use.

The MSS3 vaporiser is designed to operate with dry medical gases and is calibrated specifically for use with the anaesthetic drug range, using laser refraction technology, to ensure that all important accuracy of use.

In addition the MSS3 has safety lock function, to prevent accidental turn on - making the MSS3 not only one of the most accurately manufactured and certified vaporisers, but one of the safest.



SERVICING SPECIALIST...

With such a specialist team, dealing with a technical and refined product, MSS also provide servicing and calibration options for MSS vaporisers. In addition MSS can service and calibrate TEC 4 and TEC 5 vaporisers, plus undertake key fill/screw fill conversions - and visa versa.

At MSS, we understand the importance of accuracy and diligence in servicing and repairs undertaken. All our servicing and conversions are guaranteed and supported by MSS, giving customers peace of mind and ensuring vaporisers which leave our factory are ready to meet the demands in a critical environment.

FOR FULL TECHNICAL SPECIFICATION VISIT MSS INTERNATIONAL AT WWW.MSSMEDICAL.CO.UK

**MEETING THE
DEMANDS IN
A CRITICAL
ENVIRONMENT...**



THE COMPLETE VAPORISER PACKAGE...

As well as our world beating vaporisers, MSS International build Full Anaesthesia Machines, offering the complete package in dry medical gas delivery. Our purpose built frame enables oxygen concentrator and flow meter attachments. The MSS Anaesthesia Machines can be tailored to either a wall mounted or a mobile device.

ENGINEERED TO THE SAME STANDARDS...

In addition, MSS International also engineer and manufacture high quality parts, adapters and accessories. Designed to operate in conjunction with and alongside MSS vaporisers. For more information about parts and spares, get in touch with the team at MSS International, enquire@mssmedical.co.uk.



MSS
INTERNATIONAL LTD

**THE SPECIALIST
MANUFACTURER
OF HIGH QUALITY
ANAESTHETIC VAPORISERS**

**SERVICE AND
REPAIR**

Unit 5
Valley Road Business Park
Valley Road · Keighley
West Yorkshire · BD21 4LZ

Telephone 01535 692424
Facsimile 01535 692525

**sales@mssmedical.co.uk
www.mssmedical.co.uk**